

# ABL800 BASIC analyzer



Meeting your essential needs

# The facts about the ABL800 BASIC

## Color touch screen

- Intuitive user interface
- Stand-by option
- User-definable shortcuts
- Parameter bar provides analyzer status at a glance
- Onboard help system via multimedia tutorials

## Connectivity

- Full LIS/HIS connectivity reduces transcription errors
- Data capture through continuous synchronization
- Remote control of decentral analyzers from central location

## Inlet

- Inlet design reduces risk of air bubbles
- Automated sample aspiration
- Sample volume down to 35  $\mu\text{L}$
- High sample success rate of capillary samples

## Analytical performance

- Automated micromode reduces risk of user-induced bias and ensure accuracy of samples as small as 35  $\mu\text{L}$
- Interference-protected glucose and lactate results

## Automated quality control

- Prescheduled QC checks
- Automatic storage of patient, QC, calibration and system data
- QA Portal-compatible



## Test menu

pH |  $p\text{CO}_2$  |  $p\text{O}_2$  |  $s\text{O}_2$  | ctHb |  $\text{cK}^+$  |  $\text{cNa}^+$  |  $\text{cCa}^{2+}$  |  $\text{cCl}^-$   
cGlu | cLac

# Meeting your essential needs



The ABL800 BASIC is the ideal choice for you who are looking for efficiency and value in your critical care testing.

The ABL800 BASIC offers determination of the most important critical care parameters, including optimal values for total hemoglobin and oxygen saturation - with maximum accuracy and minimum effort.

Easy to use, the ABL800 BASIC provides you with fast and reliable results for fast patient treatment.

In addition, automatic quality control, automated sample aspiration and IT connectivity help reduce your workload and minimize errors.

# Specifications

## Default parameters

Type	Parameters	Units	Measuring ranges
pH	pH	pH scale	6.300-8.000
	cH <sup>+</sup>	nmol/L	10.0-501
Blood gas	pCO <sub>2</sub>	mmHg	5.0-250
		kPa	0.67-33.3
		Torr	5.0-250
	pO <sub>2</sub>	mmHg	0.0-800
		kPa	0.00-107
		Torr	0.0-800
Oximetry	ctHb	g/dL	0.00-27.7
		g/L	0.0-277
		mmol/L	0.00-17.2
	sO <sub>2</sub>	%	0.0-100
		Fraction	0.000-1.000

## Additional parameters

Type	Parameters	Units	Measuring ranges
Electrolyte*	cCl <sup>-</sup>	mmol/L	7-350
		meq/L	7-350
	cCa <sup>2+</sup>	mmol/L	0.20-9.99
		meq/L	0.40-19.98
		mg/dL	0.80-40.04
	cK <sup>+</sup>	mmol/L	0.5-25.0
meq/L		0.5-25.0	
cNa <sup>+</sup>	mmol/L	7-350	
	meq/L	7-350	
Metabolite*	cGlu	mmol/L	0.0-60
		mg/dL	0-1081
	cLac	mmol/L	0.0-30
		mg/dL	0-270
		meq/L	0.0-30

\* The total number of electrolytes and metabolites is limited to five.

## Measuring modes

Mode	Sample volume (µL)	Measuring time (sec)	Cycle time (sec)
All parameters	195	80	150
All parameters, micro	95	135	200
pH + BG + Oxi, micro	85	80	145
pH + BG + Oxi, micro	55	100	170
Glu + Lac, micro	35	80	145
Oxi, micro	35	80	145

Mode	Sample volume (mL)	Measuring time (sec)	Cycle time (sec)
Expired air	15	65	170

Small sample volumes.

## AutoCheck automatic quality control

Number of ampoules in carousel:	0-20
Positioning of ampoules in carousel:	Random
Lot change:	2 lots of same level possible at the same time
Liquid volume in ampoule:	700 µL
Expiration of ampoules:	24 months at 25 °C/77 °F (including 15 days at up to 32 °C/90 °F )
Conditioning time (from room temperature):	15 min with filled carousel
Scanning time:	< 30 sec with filled carousel
Cycle time:	< 2 min 40 sec
Manual QC measurement possible:	Yes
Remote control:	Remote monitoring and start of measurement, e.g. via RADIANCE

## Input parameters

Type	Definition
Patient ID	Patient identification number
Patient height	The height of the patient
Patient department	Which department the patient is from
$T$	Patient temperature
Sample type	Arterial, venous, mixed venous, capillary, prof. test, cal. verification
Patient note	Notes about the patient or sample
Patient weight	The weight of the patient
Patient Accession No.	Specific sample number
Patient age	Date of birth
Patient sex	Male or female
Draw time	When the sample was taken
Date of birth	Patient date of birth
Sample site	Not specified, brachial left/right, femoral left/right, radial left/right, finger left/right, heel left/right, umbilical cord
Patient birth weight	The weight of the newborn
Patient gestational age	Period of intrauterine fetal development from conception to birth

Type	Definition
Patient name	Name of the patient
Physician	Name of the physician
Operator	Name of the operator
Operator department	Which department the operator is from
$p50(st)$	Oxygen tension at 50 % saturation of blood at standard conditions for pH, $pCO_2$ , $FCO_{Hb}$ , $FMetHb$ , $FHbF$ at 37 °C
RQ	Respiratory quotient
$FO_2(I)$	Fraction of oxygen in dry inspired air
$\dot{Q}_t$	Cardiac output
$\dot{V}O_2$	Oxygen consumption
VCO	Volume of carbon monoxide, input value for measurement of V(B)
$sO_2(\bar{v})$	Oxygen saturation of hemoglobin in mixed venous blood
$pO_2(\bar{v})$	Oxygen tension of mixed venous blood
ctHb	Total hemoglobin concentration (if not measured)
$FCO_{Hb}(1)$	Used for determining blood volume
$FCO_{Hb}(2)$	Used for determining blood volume

Parameters can be set for mandatory input and are user-definable.



# Specifications

## Derived parameters

Type	Definition
$pH(T)$	pH of blood at patient temperature
$pCO_2(T)$	Carbon dioxide tension of blood at patient temperature
$cHCO_3^-(P)$	Concentration of hydrogen carbonate in plasma
$cBase(B)$	Concentration of titrable base of blood (actual base excess)
$cBase(B,ox)$	Actual base excess at 100 % oxygen saturation
$cBase(Ecf)$	Concentration of titrable base of extracellular fluid (standard base excess)
$cBase(Ecf,ox)$	Standard base excess at 100 % oxygen saturation
$cHCO_3^-(P,st)$	Concentration of hydrogen carbonate in plasma of standardized blood (standard bicarbonate)
$cH^+$	Concentration of hydrogen ions in blood
$cH^+(T)$	Concentration of hydrogen ions in blood at patient temperature
$ctCO_2(P)$	Concentration of total carbon dioxide in plasma
$ctCO_2(B)$	Concentration of total carbon dioxide of whole blood (CO <sub>2</sub> content)
$pH(st)$	pH of standardized blood ( $pCO_2 = 40$ mmHg)
$pO_2(T)$	Oxygen tension of blood at patient temperature
$pO_2(A)$	Oxygen tension of alveolar air
$pO_2(A,T)$	Oxygen tension of alveolar air at patient temperature
$p50$	Oxygen tension at 50 % saturation of blood
$p50(T)$	Oxygen tension at 50 % saturation of blood at patient temperature
$p50(st)$	Oxygen tension at 50 % saturation of blood at standard conditions for pH, $pCO_2$ , $FCO_{Hb}$ , $fMetHb$ , $fHbF$ at 37 °C
$pO_2(A-a)$	Difference of oxygen tension of alveolar air and arterial blood
$pO_2(A-a,T)$	Difference of oxygen tension of alveolar air and arterial blood at patient temperature
$pO_2(a/A)$	Ratio of oxygen tension of arterial blood and alveolar air
$pO_2(a/A,T)$	Ratio of oxygen tension of arterial blood and alveolar air at patient temperature

Type	Definition
$pO_2(a)/FO_2(I)$	Oxygen tension ratio of arterial blood to the fraction of oxygen in inspired air
$pO_2(a,T)/FO_2(I)$	Oxygen tension ratio of arterial blood at patient temperature to the fraction of oxygen in inspired air
$cCa^{2+}(pH=7.40)$	Concentration of ionized calcium in plasma at pH 7.40
Anion Gap( $K^+$ )	Concentration difference of $K^+ + Na^+$ and $Cl^- + HCO_3^-$
Anion Gap	Concentration difference of $Na^+$ and $Cl^- + HCO_3^-$
$\dot{D}O_2$	Oxygen delivery
Hct	Fraction of the volume of erythrocytes in the volume of whole blood
$pO_2(x)$	Oxygen extraction tension of arterial blood
$pO_2(x,T)$	Oxygen extraction tension of arterial blood at patient temperature
$ctO_2(B)$	Total oxygen concentration of blood (O <sub>2</sub> content)
$ctO_2(a-\bar{v})$	Total oxygen concentration difference between arterial and mixed venous blood
$BO_2$	Oxygen capacity of hemoglobin. The maximum concentration of oxygen bound to hemoglobin in blood, saturated so that all deoxyhemoglobin is converted to oxyhemoglobin
$ctO_2(x)$	Extractable oxygen concentration of arterial blood
$fShunt$	Volume fraction of shunted venous blood in arterial blood
$fShunt(T)$	$fShunt$ at patient temperature
RI	Respiratory Index
$RI(T)$	Respiratory Index at patient temperature
$\dot{V}O_2$	Oxygen consumption
$mOsm$	Plasma osmolarity
$Q_x$	Oxygen compensation factor of arterial blood
$\dot{Q}_t$	Cardiac output
$V(B)$	Volume of blood
$sO_2$	Saturation
$FO_2Hb$	Fraction of oxyhemoglobin in total hemoglobin in blood

## Calibration data

Automatic:	Default interval:	Interval options:
1-point cal.	4 hours	after measurement, 30 min, 1, 2, 4 hours
2-point cal.	4 hours	after measurement, 1, 2, 4 hours
System alignment	24 hours	
Cleaning	8 hours	8, 24 hours
Manual:		
tHb calibration	3 months	never, 7 days, 1, 2, 3, 4, 6 months

## IT connectivity

### Computer specifications

Intel Celeron Processor  
3 USB ports  
128 MB RAM  
Hard disk 40 GB  
TFT 10.4" VGA color touch screen  
Dedicated 80386 CPU for wet-section operations

### Data capacity

Patient results:	2000
Calibration results:	1000
QC results:	1500
System messages and service registrations:	5000

### Communication

Access to Local Area Network for backup, etc.:  
using PC network operating systems  
supporting Windows®XP

### Output protocols:

#### High-level protocols

ABL700-compatible ASTM (E1394-91)  
ABL700-compatible HL7 (Version 2.2)  
ABL5xx-compatible ASTM (E1394-91)

#### Low-level protocols

ASTM (E1381-95)  
Radiometer network protocol (TCP/IP only)  
Raw (serial only)

#### Transport layer

TCP/IP  
RS232  
RADIANCE

Communication via TCP/IP



## Simpler, faster, better

Radiometer's products and services simplify and automate all phases of acute care testing, providing you with the speed and ease of use you want and the accuracy you need.

This is acute care testing truly made simpler, faster and better.

### Contact us

Radiometer is represented in more than 100 countries.

For the local Radiometer representative nearest you, please go to [www.radiometer.com](http://www.radiometer.com).



Radiometer Medical ApS

Åkandevej 21

2700 Brønshøj

Denmark

Tel: + 45 38 27 38 27

Fax: + 45 38 27 27 27

[www.radiometer.com](http://www.radiometer.com)



Data subject to change without notice.

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, PICO and CLINITUBES are trademarks of Radiometer Medical ApS.

Windows® is a registered trademark of Microsoft Corporation.