

AQT90 FLEX analyzer



Powering your critical decisions

The AQT90 FLEX immunoassay analyzer powers your critical decisions



The AQT90 F offers a compre panel and marker infection and

Results in minutes

Short TAT facilitates timely clinical decision-making and optimizes patient flow

Ease of use

Intuitive interface and automated steps reduce demands on staff and eliminate user variability

Increased patient & operator safety

Automatic patient-sample match and no contact with blood or waste



LEX analyzer hensive cardiac s for coagulation,

d pregnancy

More control Remotely monitor and control AQT90 FLEX analyzers in real time to ensure regulatory compliance.

Better outcomes

Streamline testing while containing costs. Improve patient outcomes and satisfaction.





Ti)

High analytical performance

Results you can trust at the POC; no need for resampling



Eases remote monitoring and cost reporting

Compliance with guidelines

Compliance supported by liquid QC & QC lockout system, data capture features and comprehensive system checks

AQT90 FLEX analyzer

POWERING YOUR CRITICAL DECISIONS





ACUTE CARE TESTING

www.radiometer.com